

PRODUCT DESCRIPTION

Super **Sani-Cloth**[®] Germicidal Disposable Wipe is a premoistened nonwoven durable wipe containing a quaternary ammonium chloride/alcohol based solution. Recommended for use in hospitals and other critical care areas where the control of the hazards of cross-contamination between treated surfaces is required. Designed to be compatible with hard nonporous surfaces and equipment made of plastic, Formica[®] laminate, glass and more. Some organisms are removed from the surface by thoroughly wiping the surface with the wipe. Most remaining organisms are killed within two (2) minutes by exposure to the liquid in the wipe.

CHEMICAL COMPOSITION

Active Ingredients:

n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chlorides	0.25%
n-Alkyl (60% C_{14}^{12} , 30% C_{16}^{12} , 5% C_{12}^{12} , 5% C_{18}^{13}) dimethyl benzyl ammonium chlorides	
Isopropyl Alcohol.	
Other ingredients.	
TOTAL	
Each cloth is saturated with 5 000 parts per million of active quaternary ammonium chlorides	

Each cloth is saturated with 5,000 parts per million of active quaternary ammonium chlorides.



EFFICACY

BACTERIAL	ORGANISM	EFFICACY
MULTI-DRUG	RESISTANT B	ACTERIA:

Test Method Used: Organic Soil Load: Exposure Time: Incubation: Results:

BACTERIA:

Test Method Used: Organic Soil Load: Exposure Time: Incubation: Results:

MYCOBACTERIUM BOVIS - BCG (TB):

Test Method Used: Organic Soil Load: Exposure Time: Incubation: Results:

VIRAL ORGANISM EFFICACY ENVELOPED VIRUSES:

Test Method Used:


Organic soil load: Exposure Time: Results:

Test Method Used:

Organic soil load: Exposure Time: Results: Acinetobacter baumannii [ATCC 19606] Enterobacter cloacae NDM-1 positive [CDC 1000654] ESBL Producing Escherichia coli (E. coli) [ATCC BAA-196] Klebsiella pneumoniae KPC-2 positive, ST258 (multi-drug and carbapenem resistant) [CDC 2008030] Methicillin Resistant Staphylococcus aureus (MRSA) [ATCC 33592] Vancomycin Resistant Enterococcus faecalis (VRE) [ATCC 51299] Modified AOAC Germicidal Spray Method for Hard Surface Disinfection 5% Horse Serum or 5% Fetal Bovine Serum 2 minutes at 68-77°F 48 hours +/- 2 hours at 95-98.6°F No growth observed Bordetella pertussis [ATCC 12743] Burkholderia cepacia [ATCC 25416] Campylobacter jejuni [ATCC 29428] Escherichia coli (E. coli) [ATCC 11229] Escherichia coli 0157:H7 [ATCC 35150] Klebsiella pneumoniae [ATCC 4352] Pseudomonas aeruginosa [ATCC 15442] Salmonella enterica [ATCC 10708] Staphylococcus aureus [ATCC 6538] Modified AOAC Germicidal Spray Method for Hard Surface Disinfection 5% Horse Serum or 5% Fetal Bovine Serum 2 minutes at 66-77°F 48 hours +/- 2 hours to 6 days at 86-98.6°F No growth observed Mycobacterium bovis BCG (Tuberculosis) [ATCC 35743] Quantitative Tuberculocidal Suspension Test 5% Horse Serum 1 minute at 68°F 21 days at 98.6°F No growth observed Herpes Simplex type 2 [ATCC VR-734] Human Coronavirus [ATCC VR-740] Strain 229E Influenza A virus/Hong Kong Strain [ATCC VR-544]* * Pandemic 2009 H1N1 Influenza A virus (Kill claim included) Influenza A (H1N1) virus [ATCC VR-98] Strain A/Malaya/302/54 Vaccinia virus [ATCC VR-1354] This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the time of test for determining the virucidal efficacy of disinfectants intended for use on dry inanimate surfaces. 5% fetal bovine serum. 2 minutes at 68°F Virucidal according to the criteria established by the U.S. Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide. Respiratory Syncytial virus (RSV) This test was conducted according to U.S. Environmental Protection Agency guidelines in effect at the time of test for determining virucidal efficacy of disinfectants intended for use on dry inanimate surfaces. 5% fetal bovine serum. 2 minutes at room temperature (68°-77°F) Virucidal against Respiratory Syncytial virus (RSV) according to the criteria established by the U.S. Environmental Protection Agency for registration and labeling of a disinfectant product as a virucide.

Super Sani-Cloth

EFFICACY



Super Sani-Cloth

PATHOGENIC FUNGI EFFICACY YEAST ORGANISM:

Test Method Used: Organic Soil Load: Exposure Time: Incubation: Results: Candida albicans [ATCC 14053] Modified AOAC Germicidal Spray Method 5% Horse Serum 2 minutes at 72 · 74°F 7 days at 95 · 98.6°F No growth observed

TOXICITY

ACUTE INHALATION

Based on the inhalation test results, Super **Sani-Cloth** Germicidal Disposable Wipe has been classified as Toxicity Category IV for acute inhalation.

ACUTE ORAL TOXICITY

Based on the results of this study, Super **Sani-Cloth** Germicidal Disposable Wipe has been classified as Toxicity Category IV for acute oral toxicity.

ACUTE EYE IRRITATION

Based on the results of this study, Super **Sani-Cloth** Germicidal Disposable Wipe produced eye irritation that indicates the product would be classified as Toxicity Category II for acute eye irritation.



TOXICITY (continued)

ACUTE DERMAL TOXICITY

Based on the results of this study, Super **Sani-Cloth** Germicidal Disposable Wipe has been classified as Toxicity Category IV for dermal toxicity.

ACUTE DERMAL IRRITATION

Based on the results of primary skin irritation study, Super **Sani-Cloth** Germicidal Disposable Wipe has been classified as Toxicity Category III for dermal effects.

DERMAL SENSITIZATION

Based on the sensitization test results, Super **Sani-Cloth** Germicidal Disposable Wipe would not be considered a dermal sensitizing agent.

