

ASSEMBLY INSTRUCTIONS FOR

Bottle Hanger for Vials (#18885-18887)



Step 1. Slide the vial into the large round cutout of the bottle holder.



Step 2. Flip the small round cutout over the neck of the vial.



Step 3. Pull up on the notched portion of the large cutout.



Step 4. Hang the vial from the notched portion of the large ring on your IV Pole.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923
Web: GoHCL.com • Email: hcl@GoHCL.com