
INSTRUCTIONS FOR

Petri & Contact Plate Rack, 100mm (#21069)

Assembly:

- 1. Refer to drawing for orientation of parts.
- 2. Assemble the center divider to the base securing with one nut on each of the two threaded posts.
- 3. Assemble each of the rods to the base; he posts are keyed to fit with a specific orientation. Secure each rod with one nut on the threaded post.



Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com