

CLEANING GUIDE FOR

HCL® Label Dispensers (#12721, 12998 & 18894)

Materials:

Polycarbonate

Recommended Disinfectant Options:

20142 Sterile HCL® Presaturated Polypropylene Wipes 70% IPA, 9 x 11

17880 SatPax® Presaturated Cleanroom Wipes, Non-Sterile, 9 x 9

Alternate Disinfectant Options:

70% Isopropyl Alcohol based cleaners

3% Hydrogen Peroxide based cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.