

CLEANING GUIDE FOR

HCL® Locking Box w/ Removable Dividers (#21325)

Materials:

PETG

Recommended Disinfectant Options:

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12

19955 PREempt™ One-Step Surface Cleaner & Disinfectant, Trigger Spray, 32 oz.

Alternate Disinfectant Options:

70% Isopropyl Alcohol Based Cleaners

3% Sodium Hypochlorite Based Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.