

CLEANING GUIDE FOR

**HCL® Locking Wall Boxes (#18551, 18553, 18565, 18556, 18554, 18557, 18550,
18552, 18564, 18555, 3753, 3751-01, 19940)**

Materials:

Acrylic

Recommended Disinfectant Options:

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, Canister

19954-31 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, Case

Alternate Disinfectant Options:

3% Hydrogen Peroxide Based Cleaners

15% Sodium Hypochlorite Base Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.