

CLEANING GUIDE FOR

HCL® Stainless Steel Cart for Pyxis™ Medstation™ 4000 (#21343)

Materials:

Stainless Steel

Recommended Disinfectant Options:

20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPS, 9 x 12

19118 Sterile Isopropyl Alcohol, 16 oz.

Alternate Disinfectant Options:

70% Isopropyl Alcohol

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.