

CLEANING GUIDE FOR

HCL® Locking Refrigerator Boxes, Tilting & Non-Tilting (#18541-18544)

HCL® Three-Compartment Fridge Box Key (#18798)

Materials:

Acrylic

Recommended Disinfectant Options:

19954 PREempt™ One-Step Surface Cleaner & Disinfectant Wipes, 6 x 7, Canister

19954-31 PREempt™ One-Step Surface Cleaner & Disinfectant Wipes, 6 x 7, Canister, Case

Alternate Disinfectant Options:

3% Hydrogen Peroxide Based Cleaners

15% Sodium Hypochlorite Based Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.