

# ALCOH-Wipe<sup>®</sup> Sterile



70% USP Isopropyl Alcohol  
Saturated Clean Room Wipe  
Filtered at 0.2  $\mu$ m and  
Gamma Irradiated at 10<sup>-6</sup> SAL



## Ingredients:

USP Isopropyl Alcohol CAS# 67-63-0.....70%  
USP Water For Injection.....30%

**DANGER- FLAMMABLE**  
**KEEP OUT OF REACH OF CHILDREN**

KEEP FROM HEAT OR OPEN FLAME, DO NOT STORE ABOVE 120°F.

**FIRST AID: Eye Contact:** Immediately flush eyes with plenty of water for at least 15 minutes. Remove contact lenses, if present. Get medical attention if irritation persists after washing. **Skin Contact:** Wash skin with soap and water. Get medical attention if irritation develops and persists. **Inhalation:** Move into fresh air and keep at rest. If breathing stops, provide artificial respiration. Get medical attention if any discomfort continues.

Reorder # See Lot # Label • SDS# VEL-038 • Rev 12/13

For Spill/Exposure Emergency Response from the USA and Canada call CARECHEM 24 at 866-928-0789. From the EU call call 44 1235 239 670

**Veltek Associates, Inc.**

15 Lee Blvd.  
Malvern, PA 19355-1234 USA  
Tel: 610-644-8335  
Fax 610-644-8336  
www.sterile.com  
US Patent 6123900, 6333006, 6607695

**Veltek Associates, Inc.**

Postbus 1062  
8200 BB Lelystad, Nederland  
Tel: +800-00888700  
Fax: 610-644-8336  
www.sterile.com

Made in USA

