

# ASSEMBLY INSTRUCTIONS FOR

## Bottle Hangers for Vials



**Step 1.** Slide the vial into the large round cutout of the bottle holder.



**Step 2.** Flip the small round cutout over the neck of the vial.



**Step 3.** Pull up on the notched portion of the large cutout.



**Step 4.** Hang the vial from the notched portion of the large ring on your IV Pole.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923  
Web: GoHCL.com • Email: hcl@GoHCL.com