

Instructions for the Stainless Steel Tablet Crusher for Tall Narrow Cups (HCL #7391-01)



1. Place tablet(s) into a 1 ounce cup.



2. Place the cup into the oval area of the crusher base and cover the tablet(s) with a 2nd cup.



3. Use an up and down tapping motion to crush the tablet(s).



4. Remove the cups from the crusher. The crushed particles are ready to be dispensed.

Your small friend, indeed®
Health Care Logistics
INC.



Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923

International Call 740 477 1686 • International Fax 740 477 2923 • Web site <http://www.hcl-intl.com> • E-mail: HCL@hcl-intl.com

Instructions for the Stainless Steel Tablet Crusher for Tall Narrow Cups (HCL #7391-01)



1. Place tablet(s) into a 1 ounce cup.



2. Place the cup into the oval area of the crusher base and cover the tablet(s) with a 2nd cup.



3. Use an up and down tapping motion to crush the tablet(s).



4. Remove the cups from the crusher. The crushed particles are ready to be dispensed.

Your small friend, indeed®
Health Care Logistics
INC.



Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923

International Call 740 477 1686 • International Fax 740 477 2923 • Web site <http://www.hcl-intl.com> • E-mail: HCL@hcl-intl.com