

Instructions for #7503 Security Bag that fits Top Metro® Lifeline™ Crash Cart Box



1. Insert the crash cart box into the opening of the crash cart bag.



2. Pull off the clear plastic backing to expose the adhesive tape.



3. Seal the bag by folding over the edge of the bag and pressing it against the adhesive tape.



4. To open the bag, tear at the center perforation. Pull one half to the left and the other to the right.



HCL-INTL is a registered trademark of Health Care Logistics, Inc.

Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923
International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com