

# Instructions for #7510 Security Bag for Full Size Crash Cart Box



1. Insert the crash cart box into the opening of the crash cart bag.



2. Pull off the paper backing to expose the adhesive tape.



3. Seal the bag by folding over the edge of the bag.



4. To open the bag, remove the box from the cart drawer and tear at the center or the end perforation. Pull one half to the left and the other to the right.



HCL-INTL is a registered trademark of Health Care Logistics, Inc.

Sheet #7510  
REV.101200

**Call Free 1 800 848 1633 or 888 HCL INTL<sup>®</sup> (425 4685) • Fax Free 1 800 447 2923**  
International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: [HCL@hcl-intl.com](mailto:HCL@hcl-intl.com)