

# INSTRUCTIONS FOR

## HCL Tamper Tuf™ Vials with Caps (#7741-7749)



**Step 1.** Hold the the Tamper Tuf Vial securely, and push the thumb tab fully outward.



**Step 2.** Pull back on the tear-away strip, being careful not to pop off the top. (Hint: To start the tear pull toward the bottom of the container.)



**Step 3.** Continue pulling the tear-away strip until you reach the end.



**Step 4.** Remove the cap by grasping the ribbed edges, twisting and lifting, while holding the vial securely with your other hand.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

Web: GoHCL.com • Email: hcl@GoHCL.com

# INSTRUCTIONS FOR

## HCL Tamper Tuf™ Vials with Caps (#7741-7749)



**Step 1.** Place the vials into the holes of the holder.



**Step 2.** Start at one end of the holder and fill the first vial. Continue moving from one vial to another until all of them are filled.



**Step 3.** Slide the holder into the channel base of the LUD Press. Make sure that the tabs are facing away from the press.



**Step 4.** Using your finger, push the tab down on the holder about to be sealed.



**Step 5.** Position the tab under the press of the and push down until you hear the tab snap closed.



**Step 6.** Slide the holder forward to the next vial.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

Web: GoHCL.com • Email: hcl@GoHCL.com