

Technical Bulletin

Cleaning and Disinfecting PETG Products:

We have tested (and recommend) the following cleaners and disinfectants on PETG products.

- Sterile Alcohol 70% spray (#17979 & 18627)
- Isopropanol Alcohol 99%
- Opti-cide3 spray (#7390, 7395 & 7396)
- Infection Control Surface Wipes - individual packets and canister wipes (#7386, 7388 & 7380)
- Viraguard 70% isopropanol alcohol wipe (from Veriden Corp.)
- Dispatch Bleach disinfectant formula wipe and spray (from Caltech)

PETG Products with Imprinting:

The imprinting on the PETG was not affected by the cleaners and disinfectants we tested and recommended (under normal cleaning conditions).

Call Free 1.800.848.1633 • Fax Free 1.800.447.2923

Web: GoHCL.com • Email: hcl@GoHCL.com