CLEANING GUIDE FOR

HCL® Wipe Dispensers (#17886, 17887, 17888, 17889 & 20966)

Materials:
PETG

Recommended Disinfectant Options:
20141 Sterile HCL® Presaturated Cleanroom Wipes, 70% IPA, 9 x 12
19955 PREempt™ One-Step Surface Cleaner & Disinfectant, Tigger Spray, 32 oz.

Alternate Disinfectant Options:
70% Isopropyl Alcohol Based Cleaners
3% Sodium Hypochlorite Based Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923
GoHCL.com • hcl@GoHCL.com

REV-01
© Health Care Logistics®, Inc. 2026